
Operationsanleitung
Surgical Technique

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Achtung: Operationsanleitung beachten
Attention: Be aware of the manual of surgical technique

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Abbreviations:

AO-QC
Mini-QC
WS
as
nas
av
st
mis
-S
-T

AO-Quick Coupling
Mini-Quick Coupling
Wrench Size
angle stable
not angle stable
angle variable
self-tapping
minimal invasive
sterile
Trial

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Warning

This description does not suffice for immediate usage of the implants and instruments.
A briefing by a surgeon experienced in this field in the handling of these instruments is highly recommended.

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This document provides information about the handling of Hofer Medical Solutions implants and instruments.

This operation manual shall be considered as an addition and under no circumstances as a substitute to existing literature about surgical methods within orthopaedics and traumatology.

The content shall be regarded as a recommendation for a standardized procedure of how to apply the products without addressing the issues of any further necessary tasks, additional operative actions and possible extensions of the surgical technique.

The actual selection of the most suitable implant and its implantation method has to happen exclusively by the surgeon based on his education and the individual diagnostic findings.

All illustrations printed here have a purely symbolic character to support the description of the surgical technique and can vary.

These operation instructions don't contain any details on the use of the instruments. Corresponding documents are available in the form of


- *Usage instructions for instruments: intra- and postoperative handling*
- *Usage instructions for implants (each implant is enclosed)*

Please ensure that the diagnosis and determination of the treatment plan are left up to the surgeon.

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Hofer Medical Solutions gladly offers detailed training in safe handling and various surgical techniques.

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Tel.: +43 (0)3382 53388
Mail: office@hofer-medical.com

Please contact our 24/7 service hotline:

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The INTEOS® Small Fragment System 2.5 from Hofer Gmbh & Co KG is a multidirectional and fixed-angle Small Fragment Implant System based on the fixative internal concept.

This system generally serves to treat diaphysial, metaphsial and epiphysial fractures in small bones.

All Hofer Medical Solutions products are the result of collective development by experienced users and our own engineers. This successful collaboration leads to products that fulfil the anatomical and functional requirements of the relevant structure, from the easiest fractures to comminute fractures, on the basis of the anatomically shaped, low-profile designs as well as the various supply possibilities.

General advantages of the system:

- *Minimised irritation of ligaments and soft tissue as a result of flat plate and screw geometry and rounded edges*
- *Lowest possible dimensions of plate thickness while increasing the plate stability through the smallest possible diameter of the „screw holes“*
- *Minimally protruding screw heads*
- *High angle flexibility of up to 70° when fully locked*

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Indications

- Fractures of the proximal radius
- Osteotomies of the proximal radius

Position of Patient and Approaches

Patient Positioning:

- Standard

Approach:

- Standard

Plate Specific Details

- Perfectly anatomically pre-bent to increase the quality of the treatment and ease of use
- High rotational stability due to the exceptional plate design
- The extra flat and slim plate design allows a minimally invasive implantation, which is crucial for a quick restoration of the original mobility and function
- Should a bending nevertheless be necessary, this is possible at the discretion of the surgeon, between the screw holes on the connecting bridges
- Minimization of irritation due to low-profile design and completely retractable screw heads
- Angle stable and polyaxial locking mechanism
- Staggered screw holes in the shaft area to avoid screw collisions
- The positioning takes place near the rim of the joint surface of the radial head and underneath the annular ligament
- The plate geometry enables subchondral screw positioning for angular stable bridging of comminuted fractures
- Maximum variability for the best possible treatment, due to the compatibility with the screws and instruments of the INTEOS® 2.5 system

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Indications

- Fractures of the proximal radius
- Osteotomies of the proximal radius

Position of Patient and Approaches

Patient Positioning:

- Standard

Approach:

- Standard

Plate Specific Details

- Perfectly anatomically pre-bent to increase the quality of the treatment and ease of use
- High rotational stability due to the exceptional plate design
- The extra flat and slim plate design allows a minimally invasive implantation, which is crucial for a quick restoration of the original mobility and function
- Should a bending nevertheless be necessary, this is possible at the discretion of the surgeon, between the screw holes on the connecting bridges
- Minimization of irritation due to low-profile design and completely retractable screw heads
- Angle stable and polyaxial locking mechanism
- Staggered screw holes in the shaft area to avoid screw collisions
- The plate width is adapted to the safe zone and can also be used there
- The positioning is distal of the annular ligament and is suitable for fractures in the neck area
- The proximal screw hole arrangement allows the insertion of up to 7 locking screws to optimally support the fracture
- Maximum variability for the best possible treatment, due to the compatibility with the screws and instruments of the INTEOS® 2.5 system

Compatible to INTEOS® bone screws
Recommended screw types:

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Screw diameter 2.2 mm
Core diameter 1.6 mm
Head diameter 3 mm
Drill bit for threaded hole angle stable & polyaxial
Compatible to ForceDRIVE T6

Titan	Durchmesser	Länge	VPe
Titanium	Diameter	Length	PU
716-140-022-006	Ø 2,2 mm	6 mm	1
716-140-022-008	Ø 2,2 mm	8 mm	1
716-140-022-010	Ø 2,2 mm	10 mm	1
716-140-022-012	Ø 2,2 mm	12 mm	1
716-140-022-014	Ø 2,2 mm	14 mm	1
716-140-022-016	Ø 2,2 mm	16 mm	1
716-140-022-018	Ø 2,2 mm	18 mm	1
716-140-022-020	Ø 2,2 mm	20 mm	1
716-140-022-022	Ø 2,2 mm	22 mm	1
716-140-022-024	Ø 2,2 mm	24 mm	1
716-140-022-026	Ø 2,2 mm	26 mm	1
716-140-022-028	Ø 2,2 mm	28 mm	1
716-140-022-030	Ø 2,2 mm	30 mm	1
716-140-022-032	Ø 2,2 mm	32 mm	1

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Screw diameter 2.5 mm
Core diameter 1.9 mm
Head diameter 3.5 mm
Drill bit for threaded hole not angle stable
Compatible to ForceDRIVE T6

Titan	Durchmesser	Länge	VPe
Titanium	Diameter	Length	PU
716-145-025-020	Ø 2,5 mm	20 mm	1
716-145-025-022	Ø 2,5 mm	22 mm	1
716-145-025-024	Ø 2,5 mm	24 mm	1
716-145-025-026	Ø 2,5 mm	26 mm	1
716-145-025-028	Ø 2,5 mm	28 mm	1
716-145-025-030	Ø 2,5 mm	30 mm	1
716-145-025-032	Ø 2,5 mm	32 mm	1

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General instructions

- The implant is anatomically shaped, however, should any adjustment be necessary, such adjustments can be made with the aid of the appropriate bending pliers
- The screws in the shaft area of the bone should always be locked bicortically

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INTEOS® 2.5 Radial head plates

Two types of plates are available for the osteosynthesis of the proximal radius: the rim plate and the buttress plate.

Both plates are designed for positioning in the Hotchkiss safety zone, which is described as the area of 105° on the radius head, which is free of impingement between ulna and radius. This safety zone is located on the opposite side of tuberositas radii.

The rim plate is positioned as close as possible to the rim of the joint surface of the radius head, while the buttress plate is located distally of the annular ligament.

The radial head plates have Kirschner wire holes for the temporary Fixation with K-wires of 1.2 mm diameter.

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© 2021 Hofer GmbH & Co KG. All rights reserved.

Operation instruction booklets, handbooks, information brochures and software are protected by copyright. Copying, duplication, translation or conversion of said document through any electronic means or machine readable form, whether that be in full or in part, is not permitted without prior written authorisation from Hofer Medical Solutions.

All further rights regarding the software are defined in the provided licensing terms.

The Hofer Medical Solutions logo is a trademark of HOFER GmbH & Co KG and is registered in Austria and additional countries.

The rights to other brand and product names mentioned in this scripture are reserved by their proprietors and are accepted and recognised here. The naming of products that are not Hofer Medical Solutions products serves solely as a source of information and not of advertising. Hofer Medical Solutions accepts no liability regarding the selection, performance or efficiency of said products.

The methods and practices, etc. quoted in this scripture have been carefully tested and verified. However, should errors or suggestions of improvement be apparent, despite said verification, we politely ask you to disclose them to us. The surgical method portrayed here merely presents one possible method for the indication that is to be treated. It remains the responsibility of the relevant surgeon to either keep to the aforementioned osteosynthetic method or to vary the method according to individual needs. Direct and indirect complications as well as sequential complications are thus, in every case, excluded.

Should you have any further questions, please do not hesitate to contact us at any time.

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Ihr Hofer Medical Solutions Team

For us, a partnership is the result of a long and happy relationship in all areas of our work. Reliable bone healing for patients, optimal handling of the systems for surgeons and surgical personnel, as well as the simple preparation of the instruments, constantly strengthen this partnership.

For us, "creating" means more than just finding solutions in the form of innovative products.

It is the result of high standards, constant development, innovative products and excellent service for patients, surgeons and surgical personnel around the clock.

Please do not hesitate to contact us if you have questions about our company, our employees or our production methods.

Your Hofer Medical Solutions Team