Operationsanleitung
Surgical Technique

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Warning

*This description does not suffice for immediate usage of the implants and instruments.
A briefing by a surgeon experienced in this field in the handling of these instruments is highly recommended.*

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This document provides information about the handling of Hofer Medical Solutions implants and instruments.

This operation manual shall be considered as an addition and under no circumstances as a substitute to existing literature about surgical methods within orthopaedics and traumatology.

The content shall be regarded as a recommendation for a standardized procedure of how to apply the products without addressing the issues of any further necessary tasks, additional operative actions and possible extensions of the surgical technique.

The actual selection of the most suitable implant and its implantation method has to happen exclusively by the surgeon based on his education and the individual diagnostic findings.

All illustrations printed here have a purely symbolic character to support the description of the surgical technique and can vary.

These operation instructions don't contain any details on the use of the instruments. Corresponding documents are available in the form of


- *Usage instructions for instruments: intra- and postoperative handling*
- *Usage instructions for implants (each implant is enclosed)*

Please ensure that the diagnosis and determination of the treatment plan are left up to the surgeon.

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Hofer Medical Solutions gladly offers detailed training in safe handling and various surgical techniques.

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Tel.: +43 (0)3382 53388

Mail: office@hofer-medical.com

Please contact our 24/7 service hotline:

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The OSTYS® 4.0 cannulated screw from Hofer GmbH & Co KG (HOFER) is a system for fixing fractures to small bones and smaller bone fragments.

Cannulated screws are usually used in combination with Kirschner's wire, meaning that the exposed bone surface is hardly ever bigger than the screw diameter.

Unlike non-cannulated screws, which usually require extensive access to the bone fragments, this method is minimally invasive.

This means less stress during the operation and is therefore also a particularly gentle procedure for patients.

General benefits:

- *Self-tapping, self-drilling and screws backwards*
- *Reduced irritation of ligaments and tissue as a result of rounded screw head geometry*
- *Possibility that the whole screw head may be immersed*
- *Optimised power transmission as a result of the ForceDRIVE screw connection with "Stick-Fit" technology*

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OSTYS® 4.0 Lochschrauben KG ss sb rs T 4.0

UA-3.1 Indications

- *Fixing fractures in small bones and small bone fragments*

UA-3.2 Contraindications

Absolute Contraindications:

- *Possible sensitivity against materials used*
- *Use on the spine*

Relative contraindications:

- *Infections or inflammations (acute, chronic, local)*
- *Reduced circulation around the affected area*
- *Reduced bone stability for accurate fixation of the implant*
- *Patients whose growth plates have yet to close*
- *Patients with are unlikely to comply with postoperative rehabilitation advice*
- *Adipositas*

UA-3.3 Screw-specific details

- *Self-tapping, self-drilling and screws backwards*
- *Optimised power transmission as a result of the ForceDRIVE screw connection with "Stick-Fit" technology*
- *Excellent self-drilling function as a result of optimised rake, clearance and wedge angles, giving a more exact incision*
- *Low screw head height*

OSTYS® 4.0 Lochschrauben KG ss sb rs T 4.0

OSTYS® 4.0 Cannulated screw short thread st sd rt T 4.0

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
Screw diameter 4.0 mm
 Shank diameter 3.0 mm
 Head diameter 5.8 mm
 Cannulation 1.8 mm
 ForceDRIVE T15

Titan Titanium	Durchmesser Diameter	Länge Length	Gewindelänge Thread Length	VPe PU
717-102-040-010	Ø 4,0 mm	10 mm	5,1 mm	1
717-102-040-012	Ø 4,0 mm	12 mm	7,1 mm	1
717-102-040-014	Ø 4,0 mm	14 mm	9,1 mm	1
717-102-040-016	Ø 4,0 mm	16 mm	11,1 mm	1
717-102-040-018	Ø 4,0 mm	18 mm	13,1 mm	1
717-102-040-020	Ø 4,0 mm	20 mm	15,1 mm	1
717-102-040-022	Ø 4,0 mm	22 mm	16,0 mm	1
717-102-040-024	Ø 4,0 mm	24 mm	16,0 mm	1
717-102-040-026	Ø 4,0 mm	26 mm	16,0 mm	1
717-102-040-028	Ø 4,0 mm	28 mm	16,0 mm	1
717-102-040-030	Ø 4,0 mm	30 mm	16,0 mm	1
717-102-040-032	Ø 4,0 mm	32 mm	16,0 mm	1
717-102-040-034	Ø 4,0 mm	34 mm	16,0 mm	1
717-102-040-036	Ø 4,0 mm	36 mm	16,0 mm	1
717-102-040-038	Ø 4,0 mm	38 mm	16,0 mm	1
717-102-040-040	Ø 4,0 mm	40 mm	16,0 mm	1
717-102-040-042	Ø 4,0 mm	42 mm	16,0 mm	1
717-102-040-044	Ø 4,0 mm	44 mm	16,0 mm	1
717-102-040-046	Ø 4,0 mm	46 mm	16,0 mm	1
717-102-040-048	Ø 4,0 mm	48 mm	16,0 mm	1
717-102-040-050	Ø 4,0 mm	50 mm	16,0 mm	1
717-102-040-055	Ø 4,0 mm	55 mm	16,0 mm	1
717-102-040-060	Ø 4,0 mm	60 mm	16,0 mm	1
717-102-040-065	Ø 4,0 mm	65 mm	16,0 mm	1
717-102-040-070	Ø 4,0 mm	70 mm	16,0 mm	1

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
It is recommended to use a flat washer on severely osteoporotic bones. The contact area is increased by the washer and the risk of the screw head penetrating the osteoporotic bone.

We also offer these washers with spikes in two different sizes

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
1.) Putting on the protective sleeve:

The protective sleeve is applied to the bone. When pressed gently, the internal trocar leaves a groove on the bone which makes the introduction of the guide wires easier.

Warning: The guide wire must have a diameter of 1.6 mm or 1.5mm if the guide wire is threaded. The length must be 150 mm.

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
2.) Introducing the guide wire:

The internal trocar is removed and the guide wire is introduced with help from the protective sleeve. The depth of penetration and position of the guide wire are monitored by image converters while the guide wire is being introduced. The guide wire should represent the final position of the OSTYS cannulated screw. For additional stability in the fracture, a further guide wire can be introduced.

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
3.) Determining screw length:

After removing all of the sleeves, the OSTYS measuring scale is introduced into the drill guide. Ensure that the tip of the measuring scale is on the bone. The length of the screw can subsequently be read directly from the bottom end of the guide wire.

Warning: the reading corresponds to the actual screw length. If the screw head has not been countersunk, this must be added to the total length (3mm).

This is also the case when using washers or washers with spikes.

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
4.) Pre-drilling (optional):

If the bone is exhibiting a strong cortical layer, this can be pre-drilled with the cannulated pre-drill (Ø-2,2 mm, AO-connection).

Warning: do not drill right up to the end of the guide wire as the wire will lose its hold (stop drilling around 5 mm before the end of the wire). Furthermore, ensure that the protective sleeve is used in order to avoid irritation and damage to the surrounding tissue.

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5.) Countersink:

In order to completely countersink the cannulated screw, it is recommended to pre-drill by means of countersinking.

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6.) *Introducing the OSTYS cannulated screw:*

It is recommended to use a flat washer (e.g. with spikes) on severely osteoporotic bones. The advantage of using a flat washer is the increase in contact area, leading to a optimised distribution of compression force. Otherwise, it is recommended to countersink the bone prior to the procedure, if not using a flat washer, in order to countersink the screw and minimise the protrusion. The method without the flat washer is recommended on patients with little soft tissue coverage. Introducing the OSTYS cannulated screw is possible via guide wire.

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© 2020 Hofer GmbH & Co KG. All rights reserved.

Operation instruction booklets, handbooks, information brochures and software are protected by copyright. Copying, duplication, translation or conversion of said document through any electronic means or machine readable form, whether that be in full or in part, is not permitted without prior written authorisation from Hofer Medical Solutions.

All further rights regarding the software are defined in the provided licensing terms.

The Hofer Medical Solutions logo is a trademark of HOFER GmbH & Co KG and is registered in Austria and additional countries.

The rights to other brand and product names mentioned in this scripture are reserved by their proprietors and are accepted and recognised here. The naming of products that are not Hofer Medical Solutions products serves solely as a source of information and not of advertising. Hofer Medical Solutions accepts no liability regarding the selection, performance or efficiency of said products.

The methods and practices, etc. quoted in this scripture have been carefully tested and verified. However, should errors or suggestions of improvement be apparent, despite said verification, we politely ask you to disclose them to us. The surgical method portrayed here merely presents one possible method for the indication that is to be treated. It remains the responsibility of the relevant surgeon to either keep to the aforementioned osteosynthetic method or to vary the method according to individual needs. Direct and indirect complications as well as sequential complications are thus, in every case, excluded.

Should you have any further questions, please do not hesitate to contact us at any time.

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Ihr Hofer Medical Solutions Team

For us, a partnership is the result of a long and happy relationship in all areas of our work. Reliable bone healing for patients, optimal handling of the systems for surgeons and surgical personnel, as well as the simple preparation of the instruments, constantly strengthen this partnership.

For us, "creating" means more than just finding solutions in the form of innovative products.

It is the result of high standards, constant development, innovative products and excellent service for patients, surgeons and surgical personnel around the clock.

Please do not hesitate to contact us if you have questions about our company, our employees or our production methods.

Your Hofer Medical Solutions Team