Operationsanleitung
Surgical Technique

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Warning

*This description does not suffice for immediate usage of the implants and instruments.
A briefing by a surgeon experienced in this field in the handling of these instruments is highly recommended.*

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This document provides information about handling of Hofer Medical Solutions implants and instruments.

This operation manual is an addition and must not be considered as a substitute to existing literature about surgical methods within orthopaedics and traumatology.

The content shall be regarded as a recommendation for a standardized procedure of how to apply the products without addressing the issues of any further required measures, additional surgical interventions and possible extensions of the surgical technique.

The actual selection of the most suitable implant and its implantation method must happen exclusively by the surgeon based on his education and the individual diagnostic findings.

All illustrations printed here have a purely symbolic character to support the description of the surgical technique and can vary.

This surgical instruction doesn't contain any details on the use of the instruments. Corresponding documents are available in the form of


- *instructions for use for instruments: intra- and postoperative handling*
- *instructions for use for implants (is enclosed with each implant)*

Please note that both the diagnosis and the determination of the treatment strategy rest solely with the surgeon.

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Hofer Medical Solutions gladly offers detailed training in safe handling and various surgical techniques.

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Tel.: +43 (0)3382 53388

Mail: office@hofer-medical.com

Please contact our 24/7 service hotline:

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The OSTYS® compression screw from Hofer Medical Solutions is a system for fixing fractures to small bones and smaller bone fragments..

Through the special screw head design, irritation to soft tissue is minimised and a minimally invasive procedure is made possible..

General benefits:

- *Self-tapping, self-drilling and screws backwards*
- *Optimised power transmission as a result of the ForceDRIVE screw connection with "Stick-Fit" technology*
- *Optimised primary stability as a result of the tapered core diameter of the screw threads*
- *The head threads being in perfect coordination with the screw threads results in an optimal, prolonged compression*

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OSTYS® Compression screw kan. sb ss T 2,8

UA-3.1 Indication

- *Fixing fractures in small bones and small bone fragments*

UA-3.2 Contraindication

Absolute contraindications:

- *Possible sensitivity against materials used*
- *Use on the spine*

Relative contraindications:

- *Infections or inflammations (acute, chronic, local)*
- *Reduced circulation around the affected area*
- *Reduced bone stability for accurate fixation of the implant*
- *Inadequate soft tissue coverage*
- *Open solid fractures with inadequate possibility of cleansing*
- *Patients whose growth plates have yet to close*
- *Patients with are unlikely to comply with postoperative rehabilitation advice*
- *Obesity*

UA-3.3 Screw-specific details

- *The head threads being in perfect coordination with the screw threads results in an optimal, prolonged compression*
- *Optimised primary stability as a result of the tapered core diameter of the screw threads*
- *Excellent self-drilling function as a result of optimised rake, clearance and wedge angles, giving a more exact incision*
- *Optimised power transmission as a result of the ForceDRIVE screw connection with "Stick-Fit" technology*
- *Self-tapping, self-drilling and screws backwards*

OSTYS® Kompressionsschraube kan. sb ss T 2,8

OSTYS® Compression Screw cannulated sd st T 2,8

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
Head diameter 3.7 mm
Thread length 3.0 mm
Thread diameter 2.8 mm
Cannulation 1.1 mm
ForceDRIVE T8

Titan <i>Titanium</i>	Durchmesser <i>Diameter</i>	Länge <i>Length</i>	Gewindelänge <i>Thread Length</i>	VPe <i>PU</i>
717-110-028-012	Ø 2,8 mm	12 mm	4,6 mm	1
717-110-028-014	Ø 2,8 mm	14 mm	5,4 mm	1
717-110-028-016	Ø 2,8 mm	16 mm	6,2 mm	1
717-110-028-018	Ø 2,8 mm	18 mm	6,9 mm	1
717-110-028-020	Ø 2,8 mm	20 mm	7,7 mm	1
717-110-028-022	Ø 2,8 mm	22 mm	8,5 mm	1
717-110-028-024	Ø 2,8 mm	24 mm	9,2 mm	1
717-110-028-026	Ø 2,8 mm	26 mm	10,0 mm	1
717-110-028-028	Ø 2,8 mm	28 mm	10,8 mm	1
717-110-028-030	Ø 2,8 mm	30 mm	11,6 mm	1
717-110-028-032	Ø 2,8 mm	32 mm	12,3 mm	1
717-110-028-034	Ø 2,8 mm	34 mm	13,1 mm	1
717-110-028-036	Ø 2,8 mm	36 mm	13,9 mm	1
717-110-028-038	Ø 2,8 mm	38 mm	14,6 mm	1
717-110-028-040	Ø 2,8 mm	40 mm	15,4 mm	1

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
1.) *Introducing the guide wire:*

The guide wire is introduced with the help of the precision drill guide. The depth of penetration and position of the guide wire are monitored by image converters while the guide wire is being introduced. The guide wire should represent the final position of the OSTYS compression screw. For additional stability in the fracture, a further guide wire can be introduced.

Warning: *The guide wire must have a diameter of 1.0 mm and a length of 150 mm.*

**Navicular bone shown transparently*

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
2.) Determining screw length:

After the placement of the guide wire, the OSTYS measuring scale is introduced in to the drill guide. Ensure that the tip of the measuring scale is on the bone. The length of the screw can subsequently be read directly from the bottom end of the guide wire.

Warning: the reading corresponds to the actual screw length. If necessary, use the next shortest screw.

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
3.) Pre-drilling (Optional):

If the bone is exhibiting a strong cortical layer, this can be pre-drilled with the cannulated pre-drill (Ø-2,2 mm, AO-connection). The precision drill guide can be used with the appropriate threaded sleeve.

Warning: Do not drill right up to the end of the guide wire as the wire will lose its hold (stop drilling around 5 mm before the end of the wire). Furthermore, ensure that the protective sleeve is used in order to avoid irritation and damage to the surrounding tissue.

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4.) Introducing the OSTYS compression screw:

The "threaded sleeve" can now be removed and the compression screw attached over the protective sleeve. This must stay attached until the fracture is closed and compressed.

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© 2022 Hofer GmbH & Co KG. All rights reserved.

Operation instruction booklets, handbooks, information brochures and software are protected by copyright. Copying, duplication, translation or conversion of said document through any electronic means or machine readable form, whether that be in full or in part, is not permitted without prior written authorisation from Hofer Medical Solutions.

All further rights regarding the software are defined in the provided licensing terms.

The Hofer Medical Solutions logo is a trademark of HOFER GmbH & Co KG and is registered in Austria and additional countries.

The rights to other brand and product names mentioned in this scripture are reserved by their proprietors and are accepted and recognised here. The naming of products that are not Hofer Medical Solutions products serves solely as a source of information and not of advertising. Hofer Medical Solutions accepts no liability regarding the selection, performance or efficiency of said products.

The methods and practices, etc. quoted in this scripture have been carefully tested and verified. However, should errors or suggestions of improvement be apparent, despite said verification, we politely ask you to disclose them to us. The surgical method portrayed here merely presents one possible method for the indication that is to be treated. It remains the responsibility of the relevant surgeon to either keep to the aforementioned osteosynthetic method or to vary the method according to individual needs. Direct and indirect complications as well as sequential complications are thus, in every case, excluded.

Should you have any further questions, please do not hesitate to contact us at any time.

ÖSTERREICH

HOFER GMBH & CO KG
Jahnstrasse 10-12
A 8280 Fürstenfeld
TEL: +43(0)3382/53388
FAX: +43(0)3382/53093
office@hofer-medical.com
www.hofer-medical.com
FN21826y
UID ATU30764704

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