

HSNesin - Hofer Sliding Nail

Elastic Stable Intramedullary Nailing



Advantages of the HSN-esin Hofer Sliding Nail System

The intramedullary nailing system HSN-esin of Hofer-medical, is characterized by its particularly easy implantation technique as a result of the special shape of the runner. **Minimally invasive** implantation as well as explantation is standard.



- The runner has a broadened, rounded front end
- Optimized as well as reduced surface pressure during insertion
- Minimized risk of perforation
- Better sliding performance
- Flattened end of the sliding nail to mark the orientation of the blade
- Plugs adjusted to the different nail diameters
- Minimally invasive implantation of the plugs by means of a guided insertion device
- Fixation of plugs independent of diameter of the point of entry of the nail
- Plugs enable different protrusions and slants of the nails





Before fixation with HSN-esin



Optimized sliding performance as a result of the special design of the runner



Minimally invasive implantation of the plug with guided insertion device



Innovative plug of the HSN-Esin HOFER Sliding Nail System



AUSTRIA

HOFER GMBH & CO KG
Jahnstrasse 10-12
A 8280 Fürstenfeld

TEL: +43(0)3382/53388
FAX: +43(0)3382/53093

office@hofer-medical.com
www.hofer-medical.com

FN21826y
UID ATU30764704



GERMANY

HOFER-MEDICAL GMBH
Ruhrtalstraße 52-60
D 45239 Essen

TEL: 0800-700 80 22 (innerhalb BRD)
FAX: +49 201 56599 1069

office@hofer-medical.de
www.hofer-medical.com

HRB24645
UID DE290596924



ITALY

HOFER MEDICAL ITALIA S.R.L.
Via della Giuliana, 66
I 00195 Rome

Tel: +39 393 8950481
FAX: +39 06 81151148

office@hofer-medical.it
www.hofer-medical.it

RM-1454597
UID IT13534721009



POLAND

HOFER GMBH & CO KG SPÓŁKA
KOMANDYTOWA ODDZIAŁ W POLSCE
Radkowska 9
PL 57-402 Nowa Ruda

Tel.Kom.: +48603224112
Fax: +48 748712273

office@hofer-medical.pl
www.hofer-medical.com

REGON 141527293
UID PL1070011771

creating partnership

For us, a partnership is the result of a long and happy relationship in all areas of our work. Reliable bone healing for patients, optimal handling of the systems for surgeons and surgical personnel, as well as the simple preparation of the instruments, constantly strengthen this partnership.

For us, "creating" means more than just finding solutions in the form of innovative products.

It is the result of high standards, constant development, innovative products and excellent service for patients, surgeons and surgical personnel around the clock.

Please do not hesitate to contact us if you have questions about our company, our employees or our production methods.

Your HOFER medical solutions Team

vCard

