

## FIXUS22 Ligamentotaxis Finger Fixator



### Adjustable distraction!

Treatment of stiffness by gradual distraction of the PIP joint with the FIXUS 22 Ligamentotaxis Finger Fixator.

Dynamic distraction promotes mobility and therefore recovery.

**Give your patients the best chance of a fast and smooth recovery!**

#### ADVANTAGES

- Adjustable distraction!
- Quick and easy to use
- X-ray transparent
- Light & strong
- Patient friendly
- Single use



### ADJUSTABLE DISTRACTION!

By using threads, the distraction can be set accurately, 1 revolution of the nut is 0.5 mm. The distraction can be easily adjusted during rehabilitation.



X-Ray transparent

### QUICK AND EASY TO USE

The instrument kit makes the fixator quick and easy to install, saving time and effort.



Round shapes

### LIGHT AND STRONG

Perfect balance between light weight and strength, for optimal patient comfort.

### SINGLE USE

For a hygienic and safe treatment.



adjustable distraction

### X-RAY TRANSPARENT

The joint area remains visible on x-rays, allowing the surgeon to radiologically check the anatomy.



Instrument kit

### PATIENT FRIENDLY

Round shapes and corners without sharp parts make this product extremely comfortable and easy to wear.

### PRODUCTS

Available in 3 different sizes.

Description	Art. code
FIXUS22 Ligamentotaxis 16mm	FIX22516
FIXUS22 Ligamentotaxis 20mm	FIX22520
FIXUS22 Ligamentotaxis 22mm	FIX22522

YOUR DISTRIBUTOR: